

Product Information

T4 DNA Ligase (5U/μl)

Catalogue Number	Size
ATR-C801	100 units (20 μl reactions)
ATR-C802	200 units (40 μl reactions)

Product Description

The **T4 DNA Ligase** catalyzes the formation of a phosphodiester bond between the adjacent 5'-phosphate and 3'-hydroxyl on the blunt or cohesive end of dsDNA. It can also catalyze the linkage of RNA with ssDNA or RNA in double stranded nucleic acids. However, it cannot catalyze linkages between single stranded nucleotides. The **T4 DNA Ligase** can be used in labelling the 3'-end of RNA, cyclizing RNA and DNA oligonucleotides, cloning of cDNA, and other manipulation of nucleic acids.

Applications

- Cloning of restriction enzyme generated DNA fragments.
- Cloning of PCR products.
- Joining of double-stranded oligonucleotide linkers or adaptors to DNA.
- Site-directed mutagenesis.
- Amplified fragment length polymorphism (AFLP).
- Ligase-mediated RNA detection (3).
- Nick repair in duplex DNA, RNA or DNA/RNA hybrids.
- Self-circularization of linear DNA.

Storage

Stored at -20°C.

Shipping

The kit is shipped with ice gel transport at ≤0°C.

NOTE:

in 10 × Ligase Buffer it is normal if a small amount of precipitation occurs, please invert and mix before use.

Protocols

1. Gently vortex and briefly centrifuge **T4 DNA Ligase (5U/μl)** and **10× Ligase Buffer** after thawing.
2. Place a microcentrifuge tube on ice and add the following components for each reaction:

10 × Ligase Buffer	1 μL
Insert DNA ^a	0.3 pmol
Vector DNA ^b	0.03 pmol
T4 DNA Ligase (5U/μl)	1 μL
Water, nuclease-free	to 10 μL

a. The molar ratio of Insert/Vector should be between 3:1 and 10:1.

b. The blunt-end vector should firstly be dephosphorylated to avoid self-cycling.

3. Incubate overnight at 16°C.
4. Transformation.
5. Thaw the competent cells on ice.
6. Add the ligation product to 100 μl of competent cells. Gently flick the tube (Please do not oscillate) and incubate on ice for 30 min.
▲ The volume of the ligation product should be less than 1/10 of the volume of competent cells.
7. Incubate the mixture at 42°C in a water bath for exactly 45 sec. Then immediately chill on ice for 2 - 3 min without disturbing the mixture.
8. Add 900 μl of LB or SOC medium (without antibiotics) to the centrifuge tube. Then, shake at 37°C for 1 h at 200 - 250 rpm.
9. Warm selection plates to 37°C.
10. Centrifuge at 5,000 rpm (2,400 × g) for 5 min and discard 900 μl of supernatant. Resuspend the cells with the remaining medium and gently spread it onto a selection plate.
11. Incubate the plate upside down in a 37°C incubator for 12 - 16 h.

Precautions and Disclaimer

This product is for R&D use only, not for use in drug, diagnostic procedures, household, or other uses. When working with the product, always wear a suitable lab coat and disposable gloves, protective eyewear. When radioactive tracers are used, standard procedures for safely handling radioactive materials should be followed. For more information, please consult the appropriate material safety data sheets (MSDSs). These are available online as pdf-file or on request (info@atrmed.com). To the extent allowed by law, ATR-MED Inc. will not be liable for special, incidental, indirect, punitive, multiple, or consequential damages in connection with or arising from this document, including your use of it. By use of this product, you accept all the terms and conditions of ATR-MED products. All trademarks are the property ATR-MED unless otherwise specified.

Limited product warranty

ATR-MED® warrant its products as set forth in the ATR-MEDs' General Terms and Conditions of Sale at www.atrmed.com/terms-and-conditions.html. If you have any questions, please contact ATR-MED at www.atrmed.com/support. ATR-MED warrants that at the time of the quality release or subsequent retest date, this product confirms the information contained in this publication. The information in this guide is subject to change without notice. Purchasers must determine the suitability of the product(s) for their particular use. Please check the invoice or packing slip for additional terms and conditions for sale.

Trademarks

ATR-MED® is a registered trademark of Acell Teb Rad Med Iranian.